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PARTIAL DENTAL ATTACHMENTS 1

PRECISION FREE END PARTIAL DENTURE CONSIDERATIONS 12

RIGID PARTIAL DENTURE ATTACHMENTS 14

STERLATCH 21

BILOC AND PLASTA 23

C&M McCollum 44

RESILIENT PARTIAL DENTURE ATTACHMENTS 74

STERN ERA AND ERA-RV 76

STERN MICRO ERA 91

DALBO 104

HARDER VERTICAL 117

IC PLUNGER 124

OCTOLINK MINIATURE ANCHOR 134

CROSS-ARCH ROACH 149

SEGMENTING ATTACHMENTS 154

FIXED AND REMOVABLE BRIDGEWORK CONSIDERATIONS 158

BEYELER 158

TUBE LOCK 158

INTERLOCK 158

www.sterngold.com.cn

Alternate technique - soldering the male to the RPD framework:

1. Purchase a new Combi-Snap male and fit. Modify the male by cutting the hole up off of the locking screw. The threaded part of the screw is used to chase the screw hole in the male. With the hole tip removed the spring loads of the male can function.
2. Deactivate the male's retention using the wedge shaped end of the adjusting tool. Place the screw male in the patient's original Impression. Make an impression, picking up the male from the denture.
3. Put the rig on the male in the impression and pour the impression in dental stone. The model will have a metal representation of the Combi-Snap female.

For the soldering procedure see the technical instructions for STERN LATCH ATTACHMENTS SECTION III A, Soldering the Male to the Partial.

Note: The special Combi-Snap Screwdriver fits the slotted screw perfectly and has a plastic sleeve which can be slid over the head of the screw to hold it when the screw is being delivered to the mouth.

Combi-Snap is a registered trademark of Cendres & Minox, S.A.

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PLASTIC DOVETAIL


Summary

Bi-directional semi-precision slide attachment.
Precision retention, not adjustable.
Plastic pattern male and female.
Built in paralleling mandrels for both male and female. Allows usage with male or female positioned towards the occlusal (may be inverted for aesthetic reasons).
Can also be used in conjunction with a retentive deep set as a rigid attachment for removable partial dentures, although this is not generally recommended.

Features: Male - cast as part of crown or partial pattern.
Female - cast as part of crown pattern.

If using the attachment in a removable prosthesis, cast plastic components using alloys with a minimum Vickers hardness of 200 and at least 85,000 psi ultimate tensile strength. Appropriate choices are Pegasus ceramic alloy and Sterngold 100 crown and bridge alloy for yellow gold castings.

Minimum Space Required:				
PD	Height	FC Width	Prep Depth	RPD Height
	2.0mm	2.0mm	2.0mm	3.5mm

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4. Fabricate wax coping. If you are inexperienced with attachments, or if the abutment is tipped or rotated, it is best to create a full coverage wax pattern.

5. Place the paralleling transfer into the parallelizer. Slide the female onto the paralleling mandrel. If the wax pattern is full coverage and a box into the pattern and scribble a fine into this. This can easily be done by positioning the female beside the wax pattern and scribble a fine into the wax with a sharp instrument. Use the female as a guide. Remove enough wax from the box to prevent porcelain from contacting the attachment as deep as necessary to maintain proper tooth contour and as shallow as possible (Fig. 3).

6. Incorporate the female into the coping. If the abutment is to have a full coverage porcelain crown, the attachment must have a thin layer (4 mm minimum) of alloy surrounding the female and supporting the porcelain. If the occlusal surface of the casting is not done, checking of the porcelain can result (Fig. 4).

7. Carefully remove the paralleling transfer from the female. Place an instrument against the occlusal of the attachment as the mandrel is being withdrawn. This will help to keep the attachment in place.

8. As stated above, it is always a good idea to incorporate a lingual tracing arm into the partial frame design. The tracing arm reduces the tensile and compressive forces on the attachment, aids in lateral stability, and helps guide the main portion of the attachment into the female as the patient inserts the denture. This tracing arm becomes necessary when the attachment is reduced to 4 mm or less (Fig. 5).

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
9. Form a groove in the wax around the female face plate with a fine explorer or hand Parker (Fig. 6). Hold the Male of the faceplate, creating a bevel into which the casting investment will flow. This will help to stop alloy from casting onto the faceplate.

If you have waxed to the occlusal of the female, remove a small amount of wax from this area. This will allow the casting investment to cover the top edge of the female and prevent casting alloy from flowing inside the attachment.

10. Sprue so that the alloy will not flow directly towards the attachment. Place the pattern to provide easy access to the female when investing (Fig. 7). A direct spruing technique can be used for some single units, but an indirect technique will give you better castings.

11. The female and proximal plate must be absolutely clean to ensure an accurate casting. A cotton swab soaked in alcohol does a good job of cleaning wax and oil residue off the metal surface.

Carefully fill the female with investment, using a small instrument or brush, before the casting ring is placed onto its base (Fig. 8). Debubblers are not recommended on any part of the attachment as they may produce a film that will cause flashing.

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
STERN ERA® AND® ERA-RV

Summary

- Unilateral semi-provisional attachment.
- Nylon male, plastic pattern female.
- Female is waxed into the pattern of the abutment retention for a one-piece casting. A gingival off-set on the female allows space for the interdental papilla.
- Two offset females (2.5 and 4.5) that allow the female cycles to be placed closer to the tissue and are packaged with RV resin.
- Black lubrication male with built-in spacer.
- Six color coded males for six levels of retention: Lightest to strongest: white, orange, blue, grey, yellow, red. The white and orange are final males. The blue and grey males are over-retentive.
- Two Extra Over-retentive Males for greater retention in problem cases - yellow has a larger button than the grey, red has a larger button than the yellow.
- Nylon male greatly reduces the rate of female wear. Replacing a worn male with a new one of the same color restores the attachment's retention.

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• Optional ERA Metal jacket, or ERA-RV Metal jacket, holds the attachment male in the denture base and is sold pre-loaded with a black processing male.

• Males easily changed without use of autopolymerizing acrylic.

• Standard and reduced vertical (ERA-RV) size males. Both size males fit the same females.

Fixation: Male - retained directly in processed denture acrylic or an ERA Metal jacket. Female - cast as part of crown pattern.

Cast plastic components using alloys with a minimum Vickers hardness of 200 and at least 85,000 psi ultimate tensile strength. Appropriate choices are Pegasin ceramic alloy and Sterngold 100 green and bridge alloy for yellow gold castings.

Minimum Space Required:	Height	PC width	Prep depth	RC width
Standard	4.0mm	2.0mm	Normal crown	6.3mm
Stern ERA-RV and offsets	3.5mm	2.0mm	Normal crown	6.3mm

ERA Metal jacket thickness: 0.3mm
+ Add 1.0mm for patients with habitually strong bites.

Indications

- Unilateral or bilateral partial dentures where a resilient prosthesis is indicated.
- Used as a retentive element for overdenture bases.
- Well suited as a partial denture attachment for bonded retainers.

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